

ResQPOD[®] ITD Key Study Summary

The ResQPOD, or an earlier version of the impedance threshold device (ITD), has been the subject of over 50 published animal and clinical studies. In 2011, *The Lancet* published the first clinical trial demonstrating improved long-term survival following cardiac arrest with device technology. In this study, when the ResQPOD was used in combination with active compression decompression cardiopulmonary resuscitation (ACD-CPR), **patients had a 53% improvement in survival to hospital discharge with favorable neurologic outcome, and this survival benefit persisted to one year.**(40) An ITD carries a Class II recommendation as a CPR adjunct in the 2010 American Heart Association (AHA) guidelines.(34)

HUMAN CLINICAL TRIALS

The ResQPOD ITD has been evaluated in 18 clinical trials during both

- ◆ Conventional, standard manual CPR: 16,17,19,25,26,29,30,31,32,33,35,38,41,51
- ◆ ACD-CPR: 4,13,14,18,26,40,46,47,48

These studies have shown that the ResQPOD:

- Improves hemodynamics:
 - ◆ Increased ETCO₂ (4,38,51)
 - ◆ Systolic BP during cardiac arrest improved 20 - 97% (4,17)
 - ◆ Mean coronary perfusion pressure improved 70% (4)
- Improves short- and/or long-term survival from prehospital cardiac arrest:
 - ◆ Survival to ED admission improved 50 - 71% (19,38)
 - ◆ Survival to 24 hrs in all patients improved 45 - 68% (13,14)
 - ◆ ROSC rates improved 31 - 80% (4,25)
 - ◆ Survival to hospital discharge improved 30 - 98% (25,30,41)
 - ◆ Survival to hospital discharge with favorable neurologic outcome improved 38 - 120% (30,31,35,40,41,47) even in the absence of therapeutic hypothermia (46)
 - ◆ Six-fold improvement in survival to 90 days with favorable neurologic outcome independent of therapeutic hypothermia (48)
 - ◆ Survival to one year with favorable neurologic outcome improved 49% (40)
 - ◆ Meta-analysis showed more than doubling of favorable neurologic outcome (26)
- Improves short- and/or long-term survival from in-hospital cardiac arrest:
 - ◆ Survival to hospital discharge rates improved 60 - 65% with adoption of AHA guidelines (including an ITD) (32,33)
- Provides benefit in non-V-fib cardiac arrest rhythms:
 - ◆ In PEA patients, survival to 24 hrs more than doubled (16) and survival to hospital discharge improved >100% (32)
 - ◆ Survival in patients presenting in asystole tripled (19)
- Works effectively on a variety of airway adjuncts (3,18,36 [manikin])
- Is clinically and cost-effective (41)

ANIMAL STUDIES

The ResQPOD ITD has been evaluated in 23 animal studies during both

- ◆ Conventional, standard manual CPR: 2,3,5,7,8,12,15,22,23,24,25,27,50
- ◆ ACD-CPR: 1,3,6,7,9,10,11,12,20,21,24,28,37,39,42,43,44,45,49,50,52,53

These studies have shown that the ResQPOD:

- Improves hemodynamics and vital organ blood flow:
 - ◆ Increased cardiac output (22), coronary perfusion pressure (1,5,9,10,15,20,21,22,24,25,50) and blood flow to the heart (2,5,7[doubles],10)
 - ◆ Increased cerebral perfusion pressure (20,22,23,24,25,49,50) and blood flow to the brain (1,2[≥50%],5,9,10,20,22,23,24,25,27,49)
 - ◆ Raised aortic blood pressure (1,8,9,10,11,15,24,50)
 - ◆ Lowered intracranial pressure during the chest wall recoil phase of CPR (23,25,49)
 - ◆ Increased ETCO₂ (8,22)
- Enhances negative intrathoracic pressure with an LMA (3,12)
- Improves survival (1,8,21,28) and neurologically-intact survival (8,20,49,53)
- Despite prolonged downtimes, improves hemodynamics and survival when used in combination with sodium nitroprusside (37,39,44,45) with good neurologic recovery (39,42,43,52)
- Improves cerebral metabolism (11) and hemodynamics (9,11) during hypothermic cardiac arrest, and induces cerebral hypothermia more rapidly after ROSC (21)
- Increases the likelihood of successful defibrillation (1) or the total energy required for successful defibrillation (21,23,28)
- Circulates drugs more effectively (9)
- Improves hemodynamics in a pediatric model of cardiac arrest (10,23,25)
- Optimizes and complements current AHA CPR recommendations (22,27)

Finally, the best outcomes following cardiac arrest will be achieved combining a continuum of care and therapies, not a single drug or device. ACSI supports the approach taken by the Take Heart America™ Demonstration Project, which promotes a full spectrum of optimal therapies, including public recognition, widespread CPR training, performance of high-quality CPR with an ITD, and definitive, specialized care at Level One Cardiac Arrest Centers (41) offering state-of-the-art post-resuscitation care to optimize neurologic recovery (e.g. therapeutic hypothermia). Go to www.takeheartamerica.org for more information.

IMPEDANCE THRESHOLD DEVICE BIBLIOGRAPHY

1. **Lurie KG** et al. Improving ACD-CPR with an inspiratory impedance valve. *Circulation* 1995;91(6):1629-1632.
2. **Lurie KG** et al. Optimizing standard CPR with an inspiratory impedance threshold valve. *Chest* 1998;113(4):1084-1090.
3. **Lurie KG** et al. Potential role of vocal cords during CPR. *Circulation* 1998;98(17):1478.
4. **Plaisance P** et al. Inspiratory impedance during ACD-CPR: a randomized evaluation in patients in cardiac arrest. *Circulation* 2000;101(9):989-994.
5. **Lurie KG** et al. Improving standard CPR with an inspiratory impedance threshold valve in a porcine model of cardiac arrest. *Anesth Analg* 2001;93(3):649-655.
6. **Voelckel WG** et al. The effects of PEEP during ACD-CPR with the inspiratory threshold valve. *Anesth Analg* 2001;92(4):967-974.
7. **Langhelle A** et al. Inspiratory impedance threshold valve during CPR. *Resuscitation* 2002;52:39-48.
8. **Lurie KG** et al. Use of an inspiratory impedance valve improves neurologically intact survival in a porcine model of ventricular fibrillation. *Circulation* 2002;105(1):124-129.
9. **Raedler C** et al. Vasopressor response in a porcine model of hypothermic cardiac arrest improved with ACD-CPR using the inspiratory impedance threshold valve. *Anesth Analg* 2002;95(6):1496-1502.
10. **Voelckel WG** et al. Effects of ACD-CPR with the inspiratory threshold valve in a young porcine model of cardiac arrest. *Pediatr Res* 2002;51(4):523-527.
11. **Bahlmann L** et al. Brain metabolism during CPR assessed with microdialysis. *Resuscitation* 2003;59(2):255-260.
12. **Lurie KG** et al. Evaluation of a prototypic inspiratory impedance threshold valve designed to enhance the efficiency of CPR. *Respir Care* 2003;48(1):52-57.
13. **Wolcke BB** et al. Comparison of standard CPR versus the combination of ACD-CPR and an inspiratory impedance threshold device for out-of-hospital cardiac arrest. *Circulation* 2003;108(18):2201-2205.
14. **Plaisance P** et al. Evaluation of an ITD in patients receiving ACD-CPR for out of hospital cardiac arrest. *Resuscitation* 2004;61:265-271.
15. **Yannopoulos D** et al. Reducing ventilation frequency combined with an inspiratory impedance device improves CPR efficiency in a swine model of cardiac arrest. *Resuscitation* 2004 ;61:75-82.
16. **Aufderheide TP** et al. Clinical evaluation of an inspiratory impedance threshold device during standard CPR in patients with out of hospital cardiac arrest. *Crit Care Med* 2005;33:734-740.
17. **Pirralo RG** et al. Effect of an inspiratory impedance threshold device on hemodynamics during conventional manual CPR. *Resuscitation* 2005;66:13-20.
18. **Plaisance P** et al. Use of an inspiratory ITD on a facemask and ET tube to reduce intrathoracic pressures during the decompression phase of ACD-CPR. *Crit Care Med* 2005;33(5):990-994.
19. **Thayne R** et al. Use of an ITD improves short-term outcomes following out-of-hospital cardiac arrest. *Resuscitation* 2005;67(1):103-108.
20. **Metzger AK** et al. Effect of an ITD and a novel ACD-CPR device on cerebral perfusion pressures and 24-hour neurological survival in a porcine model of cardiac arrest. *Circulation* 2006;114(18):II-554.
21. **Srinivasan V** et al. Rapid induction of cerebral hypothermia is enhanced with active compression decompression plus inspiratory ITD CPR in a porcine model of cardiac arrest. *J Am Coll Cardiol* 2006;47(4):835-841.
22. **Yannopoulos D** et al. Clinical and hemodynamic comparison of 15:2 and 30:2 compression-to-ventilation ratios for CPR. *Crit Care Med* 2006;34(5):1444-1449.
23. **Alexander C** et al. Dual mechanism of blood flow augmentation to the brain using an ITD in a pediatric model of cardiac arrest. *Circulation* 2007;116(16):II-433.
24. **Yannopoulos D** et al. Lower extremity counterpulsation during the decompression phase of CPR improves hemodynamics and provides continuous forward carotid blood flow. *Circulation* 2007;116(16):II-485.
25. **Aufderheide TP** et al. From laboratory science to six emergency medical services systems: new understanding of the physiology of CPR increases survival rates after cardiac arrest. *Crit Care Med* 2008;36:S397-S404.
26. **Cabrini L** et al. Impact of ITDs on CPR: a systematic review and meta-analysis of randomized controlled studies. *Crit Care Med* 2008;36(5):1625-1632.
27. **Lurie KG** et al. Comparison of 10-breaths-per-min versus a 2-breaths-per-min strategy during CPR in a porcine model of cardiac arrest. *Respir Care* 2008;52(7):862-870.
28. **Matsuura T** et al. An ITD combined with an automated ACD-CPR device (LUCAS) improves chances for survival in pigs in cardiac arrest. *Circulation* 2008;118:S1449-50.
29. **Aufderheide TP** et al. Ventilation rate and use of the ITD are correlated with hemodynamics during CPR in humans. *Circulation* 2009;120:S-669.
30. **Aufderheide TP** et al. Implementing the 2005 AHA guidelines improves outcomes after out-of-hospital cardiac arrest. *Heart Rhythm* 2010;9(10):1357-1364.
31. **Hinchey PR** et al. Improved out-of-hospital cardiac arrest survival after the sequential implementation of the 2005 AHA guidelines for compressions, ventilations and induced hypothermia: the Wake County experience. *Ann Emerg Med* 2010;56(4):358-361.
32. **Thigpen K** et al. Implementing the 2005 AHA guidelines, including use of an ITD, improves hospital discharge rates after in-hospital cardiac arrest. *Respir Care* 2010;55(8):1014-1019.
33. **Thigpen K** et al. Implementation of the 2005 AHA guidelines improves in-hospital cardiac arrest survival rates in a community hospital: a 5-year case series. *Circulation* 2010;122:A48.
34. **2010 American Heart Association guidelines for CPR and ECC science.** *Circulation* 2010;122:S722.
35. **Dailey M** et al. Implementation of the AHA guidelines with a systems-based approach improves survival to hospital discharge following cardiac arrest. *Circulation* 2010;122:A51.
36. **Genzwerker HV** et al. Influence of an impedance threshold valve on ventilation with supraglottic airway device during CPR in a manikin. *Resuscitation* 2010;81:1010-1013.
37. **Yannopoulos D** et al. Effects of epinephrine and sodium nitroprusside on left ventricular wall thickness and cavity size and carotid blood flow during CPR. *Circulation* 2010;122:A72.
38. **Saussey J** et al. Optimization of CPR with an ITD, automated compression CPR and post-resuscitation in-the-field hypothermia improves short-term outcomes following cardiac arrest. *Circulation* 2010;122:A256.
39. **Yannopoulos D** et al. Sodium nitroprusside enhanced CPR improves survival with good neurologic function in a porcine model of prolonged cardiac arrest. *Crit Care Med* 2011;39(6):1269-1274.
40. **Aufderheide T** et al. Standard CPR versus ACD-CPR with augmentation of negative intrathoracic pressure for out-of-hospital cardiac arrest: a randomised trial. *Lancet* 2011;377(9762):301-311.
41. **Lick CJ** et al. Take Heart America™: a comprehensive, community-wide, systems-based approach to the treatment of cardiac arrest. *Crit Care Med* 2011;39(1):26-33.
42. **Yannopoulos D** et al. Sodium nitroprusside improves carotid blood flow and 24-hour survival in a porcine model of prolonged CPR. *Circulation* 2010;122:A164.
43. **Schultz J** et al. Sodium nitroprusside enhanced CPR prevents post-resuscitation left ventricular dysfunction and improves 24-hour survival and neurological function in a porcine model of prolonged untreated V-fib. *Resuscitation* 2011;82S:S35-S40.
44. **Schultz J** et al. Sodium nitroprusside-enhanced CPR improves resuscitation rates after prolonged untreated cardiac arrest in two porcine models. *Crit Care Med* 2011;39(12):2705-2710.
45. **Schultz J** et al. Sodium nitroprusside enhanced CPR (SNPeCPR) improves vital organ perfusion pressures and carotid blood flow in a porcine model of cardiac arrest. *Resuscitation* 2011;Aug 22 Epub ahead of print.
46. **Frascone RJ** et al. ACD-CPR and augmentation of negative intrathoracic pressure is neuro-protective in patients with out-of-hospital cardiac arrest. *Resuscitation* 2011;82(Suppl 1):S3.
47. **Mahoney B** et al. Use of an ITD with ACD-CPR improves survival with good neurologic function following cardiac arrest from non-traumatic etiologies. *Resuscitation* 2011; 82(Suppl 1):S1-S2.
48. **Wayne M** et al. Improvement of long-term neurological function after sudden cardiac death and resuscitation: impact of CPR method and post-resuscitation care. *Resuscitation* 2011; 82(Suppl 1):S3.
49. **Metzger A** et al. Improved cerebral perfusion pressures and 24-hour neurological survival in a porcine model of cardiac arrest with ACD-CPR and augmentation of negative intrathoracic pressure. *Crit Care Med* 2012; in press.
50. **Burnett A** et al. Potential negative effects of epinephrine on carotid blood flow and ETCO2 during ACD-CPR utilizing an ITD. *Circulation* 2011;124:A189.
51. **Chase D** et al. Effect of early use of a King airway and ITD on circulation, survival and neurologic outcome in out-of-hospital cardiac arrest. *Prehosp Emerg Care* 2012;16(1):178.
52. **Yannopoulos D** et al. Controlled pauses at the initiation of sodium nitroprusside-enhanced CPR facilitate neurological and cardiac recovery after 15 minutes of untreated ventricular fibrillation. *Crit Care Med* 2012;40(5);in press.
53. **Pantazopoulos I** et al. Use of the impedance threshold device improves survival rate and neurological outcome in a swine model of asphyxial cardiac arrest. *Crit Care Med* 2012;40(2);in press.

The generally cleared indication for the ResQPOD available for sale in the United States (US) is for a temporary increase in blood circulation during emergency care, hospital, clinic and home use. The version of the ResQPOD used in the ROC PRIMED Study and ResQ Trial is not yet approved for sale in the US. Research is ongoing in the US to evaluate the long-term benefit of the ResQPOD for indications related to patients suffering from cardiac arrest. The studies listed here are not intended to imply specific outcome-based claims not yet cleared by the US FDA.

ADVANCED CIRCULATORY SYSTEMS, INC

www.advancedcirculatory.com

1-877-737-7763

Perfusion on Demand



Mar 2012; 49-0341-000, 08